v i v a n t i s

Nucleotides

Product Datasheet

10mM dNTP Mix Product No: NP2409 Quantity: 0.25ml

Lot : Expiry date :

Concentration: 10mM

Store at -20°C

info@vivantechnologies.com

Description:
dNTP Mix is an aqueous solution containing

dATP, dCTP, dGTP and dTTP, each in a final concentration of 10mM.

Quality Control:

Functionally tested in PCR with *Taq* and *Pfu* DNA Polymerases. Purity of each dNTP >98% by HPLC.



Product Use Limitation

This product is for research purposes and in vitro use only.

V i V a n t i S | www.vivantechnologies.com

SNP2409 rev0 010311

v i v a n t i s

Nucleotides

Product Datasheet

10mM dNTP Mix

Product No: NP2410 Quantity: 1ml

Lot : Expiry date :

Concentration: 10mM

Store at -20°C

info@vivantechnologies.com

Description:dNTP Mix is an aqueous solution containing dATP, dCTP, dGTP and dTTP, each in a final concentration of 10mM.

Quality Control:

Functionally tested in PCR with *Taq* and *Pfu* DNA Polymerases. Purity of each dNTP >98% by HPLC.



Product Use Limitation

This product is for research purposes and in vitro use only.

V i V a n t i S | www.vivantechnologies.com

v i v a n t i s

Nucleotides

Product Datasheet

2mM dNTP Product No: NP2411 Quantity: 1ml

Lot : Expiry date : Concentration : 2mM

Store at -20°C info@vivantechnologies.com Description:

dNTP Mix is an aqueous solution containing dATP, dCTP, dGTP and dTTP, each in a final concentration of 2mM.

Quality Control:

Functionally tested in PCR with $\it Taq$ and $\it Pfu$ DNA Polymerases. Purity of each dNTP >98% by HPLC.



Product Use Limitation

This product is for research purposes and in vitro use only.

V i V a n t i S | www.vivantechnologies.com

DSNP2411_rev0_010311

v i v a n t i s

Nucleotides

Product Datasheet

2mM dNTP Product No : NP2412 Quantity : 5 x 1ml

Lot :
Expiry date :
Concentration : 2mM

Store at -20°C

info@vivantechnologies.com

Description:

dNTP Mix is an aqueous solution containing dATP, dCTP, dGTP and dTTP, each in a final concentration of 2mM.

Quality Control:

Functionally tested in PCR with *Taq* and *Pfu* DNA Polymerases. Purity of each dNTP >98% by HPLC.



Product Use Limitation

This product is for research purposes and in vitro use only.

V I V a 11 t I S | www.vivantechnologies.com

v i v a n t i s

Nucleotides

Product Datasheet

25mM dNTP Mix Product No: NP2413 Quantity: 0.25ml

Lot : Expiry date :

Concentration: 25mM

Store at -20°C

info@vivantechnologies.com

Description:

dNTP Mix is an aqueous solution containing dATP, dCTP, dGTP and dTTP, each in a final concentration of 25mM.

Quality Control:

Functionally tested in PCR with *Taq* and *Pfu* DNA Polymerases. Purity of each dNTP >98% by HPLC.



Product Use Limitation

This product is for research purposes and in vitro use only.

V I V a n t I S | www.vivantechnologies.com

v i v a n t i s

Nucleotides

Product Datasheet

25mM dNTP Mix

Product No : NP2414 Quantity : 5 x 0.25ml

Lot : Expiry date :

Concentration: 25mM

Store at -20°C Description:

info@vivantechnologies.com

dNTP Mix is an aqueous solution containing dATP, dCTP, dGTP and dTTP, each in a final concentration of 25mM.

Quality Control:

Functionally tested in PCR with *Taq* and *Pfu* DNA Polymerases. Purity of each dNTP >98% by HPLC.



Product Use Limitation

This product is for research purposes and in vitro use only.

V V ant V sww.vivantechnologies.com